

1/16




FIG. 1




FIG. 2




FIG. 3A




FIG. 3B

4/16




FIG. 3C




FIG. 3E










FIG. 3D

5/16



6/16



7/16




FIG. 6

8/16

FIG. 7





FIG. 8A

9/16

RECEIVER

40

FIG. 8B





FIG. 9



10/16

CORE Configuration

Module Settings

PVT Path:

Admin Capcode:

Keep Events For: (1-180 days)

Default Scheduled Sample Duration: (0-480 Minutes)

Reminder Control

- ☒ Allow Reminders for Critical Alarms
- ☐ Allow Reminders for Warning Alarms
- ☐ Allow Reminders for Advisory Alarms
- ☐ Send Reminders Continuously
- ☒ Send Scheduled Samples
- ☐ Send Highest Priority Alarms to Admin Capcode
- ☒ Send Workstation Notification

Device State Time Out: (1-30 Minutes)

Guard Time Duration: Milliseconds

Protected Period: (0-60 Seconds)

Acknowledgment Period: (15-120 Seconds)

OK Cancel Apply

Callout numbers: 910 points to the Reminder Control section; 912 points to the PVT Path field; 914 points to the Keep Events For field; 916 points to the Send Scheduled Samples checkbox; 918 points to the Send Workstation Notification checkbox; 920 points to the Protected Period field; 922 points to the Acknowledgment Period field.

FIG. 10A

FIG. 10B

11/16

The SDN Configuration window has tabs for Module, Settings, Mode, Bed Selection, Yellow Alarms, and Inop Alarms. The Bed Selection tab is active. It contains two sections for SDN Card #1 and SDN Card #2. Each section has a 'Select All' and 'Clear All' button. Below these are 24 checkboxes labeled Bed #1 through Bed #24. In the SDN Card #1 section, Bed #17 is checked. In the SDN Card #2 section, all checkboxes are unchecked. At the bottom are OK, Cancel, and Apply buttons.

944 → ☒ SDN Card #1 Select All Clear All

940 → ☐ Bed #1 ☐ Bed #7 ☐ Bed #13 ☐ Bed #19
☐ Bed #2 ☐ Bed #8 ☐ Bed #14 ☐ Bed #20
☐ Bed #3 ☐ Bed #9 ☐ Bed #15 ☐ Bed #21
☐ Bed #4 ☐ Bed #10 ☐ Bed #16 ☐ Bed #22
☐ Bed #5 ☐ Bed #11 ☒ Bed #17 ☐ Bed #23
☐ Bed #6 ☐ Bed #12 ☐ Bed #18 ☐ Bed #24

942 → ☐ SDN Card #2 Select All Clear All

941 → ☐ Bed #1 ☐ Bed #7 ☐ Bed #13 ☐ Bed #19
☐ Bed #2 ☐ Bed #8 ☐ Bed #14 ☐ Bed #20
☐ Bed #3 ☐ Bed #9 ☐ Bed #15 ☐ Bed #21
☐ Bed #4 ☐ Bed #10 ☐ Bed #16 ☐ Bed #22
☐ Bed #5 ☐ Bed #11 ☐ Bed #17 ☐ Bed #23
☐ Bed #6 ☐ Bed #12 ☐ Bed #18 ☐ Bed #24

943 →

OK Cancel Apply

FIG. 10C

The SDN Configuration window has tabs for Module, Devices, and Settings. The Devices tab is active. It contains a table with columns Name, Port, Unit, and Rank. Below the table are Add, Modify, and Delete buttons. At the bottom are OK, Cancel, and Apply buttons.

970 →

974 →

Name	Port	Unit	Rank
<input checked="" type="checkbox"/> New Device	2		2
<input type="checkbox"/> Care Unit '8	-		1
<input type="checkbox"/> Cell Network	-		4
<input type="checkbox"/> Hospital Network	-		3

960 →

968 →

972 →

969 →

962 → Add Modify Delete

964 →

966 →

OK Cancel Apply

FIG. 11A

12/16

1414

1402

1404

1406

1408

1410

1412

1416

1418

PATIENTS CAREGIVERS UNIT: LAB1 Caregiver Administration

[Caregiver Setup (Add/Remove)] [Clear Assignments For All Caregivers] [Group Setup]

Name	Receiver Number	Unit	Pendant Id	Options			
Caregiver_3	0020126	LAB1	eeddccc	[Send Page]	[Change Name]	[Patient Assignment]	[Schedules]
Caregiver_2	0020125	LAB1	aabbccc	[Send Page]	[Change Name]	[Patient Assignment]	[Schedules]
Caregiver_1	0020123	LAB1	abcdeff	[Send Page]	[Change Name]	[Patient Assignment]	[Schedules]
Everyone	0000009	LAB1		[Send Page]			
Group_1		LAB1		[Send Page]	[Change Name]	[Patient Assignment]	

FIG. 11B

1420

1404

1406

1408

1410

1412

1421

PATIENTS CAREGIVERS Caregiver Setup

[Add Caregiver]

Name	Receiver Number	Unit	Pendant Id	Options
Caregiver_1	0020123	LAB1	abcdeff	[Remove] [Edit]
Caregiver_2	0020125	LAB1	aabbccc	[Remove] [Edit]
Caregiver_3	0020126	LAB1	eeddccc	[Remove] [Edit]

FIG. 11C

1424

1426

1428

1430

1427

1432

1434

1435

PATIENTS CAREGIVERS UNIT: LAB1 Add Caregiver

Name: Jane Wilson (Required. Can be a maximum of 20 characters)

Receiver Id: 1234567

Unit: LAB1

Pendant Id: aabbccc (Optional. Must be 5 characters)

[Add] [Reset]

[Back to Caregiver Setup] [Back to Caregiver Administration]

[Set Info for Caregiver]

Groups

Available Groups

LAB1 All

Nurse Manager LAB1

Nurse Manager All

Cardiac Care

Assign

FIG. 11D

13/16

PATIENTS CAREGIVERS UNIT: LAB8 Assign Patient(s) to Caregiver_3

1438

Patients currently assigned to Caregiver_3

Position	Patient	Bed	Options
This caregiver is not assigned to any patients			

[Remove All Assignments for Caregiver_3]

Available Patients:

[Sort by Name] [Sort by Bed]

1442

Bed: BED68 (LAB8; BED68)
Bed: BED2 (Patient 2)
Bed: BED4 (8 Patient)
Bed: BED5 (Patient 7)

1440

Assign Get

1446

FIG. 11E

PATIENTS CAREGIVERS UNIT: LAB8 Assign Patient(s) to Caregiver_3

1438

Patients currently assigned to Caregiver_3

Position	Patient	Bed	Options
1	PATIENT 9	BED2	Remove Assignment
1	8 PATIENT	BED4	Remove Assignment
2 (Cardiac)	PATIENT 7	BED5	Remove Assignment

[Remove All Assignments for Caregiver_3]

1456

1450 1452 1448 1454

Available Patients:

[Sort by Name] [Sort by Bed]

1442

Bed: BED68 (LAB8; BED68)

1440

Assign Get

1446

FIG. 11F

PATIENTS CAREGIVERS UNIT: LAB8 Assign a Caregiver to this Patient

1470

Patient	Unit	Bed
Patient 8	LAB8	BED2

1474

Name	Receiver Number	Unit	Options
Caregiver_99	1234567	LAB8	Assign (1) (2) (3)
Caregiver_1	0020123	LAB8	Assign (1) (2) (3)
Caregiver_2	0020125	LAB8	Assign (1) (2) (3)
Caregiver_3	0020126	LAB8	Assign (1) (2) (3)
LAB8 All		LAB8	Assign (1) (2) (3)
Nurse Manager All		Multiple	Assign (1) (2) (3)

1472

14/16

MEI Configuration

Module Alarm Filters Settings **Parameter Alarms**

Place a check mark next to the alarms that should be transmitted as warnings.

Alarm	Crisis	Warning
ECG	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>
ST Segment	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>
NBP	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>
SpO2	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>
ART	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>
RESP	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>

OK Cancel Apply

FIG. 12B

MEI Configuration

Module Alarm Filters **Settings** Parameter Alarms

Checked Alarms are Processed.

Critical ECG Alarms will be Processed Regardless of Entries Below.

Asterisks (*) can be used for Wild Card Matching.

Add New Filter... **Delete Selected Filters**

Filtered Alarm

- ☒ ACC VENT
- ☒ APNEA
- ☒ ARR OFF
- ☒ ARR SUSPEND
- ☒ ARTIFACT
- ☒ BIGEMINY
- ☒ BRADY
- ☒ COUPLET
- ☒ HR HI*
- ☒ HR LO*
- ☒ IRREGULAR
- ☒ LEADS FAIL
- ☒ NBP HI*
- ☒ NBP LO*
- ☒ NO TELEM
- ☒ NOISE
- ☒ OFF NETWORK

OK Cancel Apply

FIG. 12A

15/16

FIG. 13




FIG. 14





FIG. 15

16/16

FIG. 16





FIG. 17